

# Physics Tutorial 3: Manipulating Contrast in MRI

*Originally Written by Dan Gallichan*

*Last Modified: Eleanor Berry, Mark Chiew, Hannah Hare, Tom Okell, Oct 2014*

The purpose of this tutorial is to attempt to provide a basic understanding of what influences some of the different pulse-sequence parameters – primarily TR (repetition time), TE (echo time) and flip angle – have on the resulting images. Questions without any asterisks are those that should be attempted by everyone, whereas more challenging questions are marked with one (\*) or two (\*\*) asterisks, and should be considered optional. Take these opportunities to think about these questions, and then discuss your answers with your tutor.

At the end of the tutorial period, you will receive a “take-home tutor”, which is an annotated version of this tutorial guide that will help you complete the tutorial at home if you don’t manage to make it all the way through with your tutor, or if you missed the tutorial session.

Attendance will be taken by the tutors, and marks will be given by participation. If you would like additional feedback or clarification on the tutorial material, you are welcome to submit your questions or comments to Weblearn, and a tutor will provide written feedback for you. Those unable to attend the tutorial must submit answers to all unstarred questions to receive credit for the tutorial.

---

## Part 0 – Getting Started

### 0.1 Starting MATLAB

Download and unzip the file containing the tutorial resources from Weblearn, or if you have access to the FMRI internal network, copy them into your current directory from here:

~mchiew/GradCourse/3\_Contrast\_Manipulation

Start MATLAB, and make sure you’re inside the tutorial directory (i.e., the folder containing all the tutorial files).

Note to jalapeno00 users: please start with the -nojvm option, “matlab -nojvm”; *this should reduce server CPU load if lots of people are trying to use jalapeno00 simultaneously*

---

## Part 1 – T1 Contrast [40 minutes]

Approximate Relaxation Times at 1.5T		
	T1 (ms)	T2 (ms)
White Matter (WM)	600	80
Grey Matter (GM)	900	100
Cerebrospinal Fluid (CSF)	3500	2000

Table 1

In this tutorial, we will be making extensive use of `simContrast`, which is a

simplified simulator for the formation of image contrast in MRI. By default, it assumes approximate values for the different tissue relaxation rates at 1.5 T, as shown in Table 1.

To use the simulator, you must specify a TR, TE and flip angle:

```
>> TR = 8000;  
>> flipAngle = 90;  
>> TE = 5;  
>> simContrast
```

*Discuss the inputs, and the various components of the simContrast figure window.*

Remember that T1 (sometimes also seen as  $T_1$ ) is a property specific to each tissue type that dictates how quickly magnetization recovers along the z (i.e. main magnetic field  $B_0$ ) direction. Tissues with larger (longer) T1 relaxation times take longer to return to equilibrium than tissues with smaller (shorter) T1's.

#### **Extra Information**

T1 recovery is exponential, and that T1 is the time constant of that exponential recovery. For example,  $M_z = M_0(1 - e^{-t/T_1})$  following a 90 degree excitation. That is, after a time T1, the magnetisation  $M_z$  would have recovered  $100 \cdot (1 - e^{-1})\%$  or ~63% of its equilibrium value.

*Demonstrate this difference by simulating the recovery curves using a long TR (e.g. TR=15000) and a large flip angle (e.g. flipAngle=90). Use any short TE  $\geq 0$ . Describe the resulting subplot of  $M_z$  in the top left of the figure window, taking care to link the shorter recovery times of the WM and GM (blue and green curves) with the small T1 values in Table 1, and the long recovery time of the CSF (red curve) with its large T1.*

Recall that the quantity  $M_z$  represents is the amount of potential or stored magnetization that is available for subsequent measurements. In terms of the battery analogy, batteries with larger T1 values take longer to recharge (recover) than batteries with shorter T1's. Therefore, the shorter the T1, the shorter the duration needed (TR) for the magnetisation to recharge in order to achieve the maximum possible signal.

T1-weighting refers to the act of producing images that reflect differences in tissue T1. To do this, the TR is often selected to maximise or accentuate the difference in signal levels associated with respective T1 values.

*Encourage students to spend a couple minutes playing with the simContrast simulator modifying TR values from extremely short to extremely long, and moderate values in between to explore the effect on measured signal strengths. Instruct them to keep flipAngle=90, and TE = 0 to simplify interpretation.*

#### **Question 1.1**

Images produced with very long TRs (and short TEs), have maximum signal levels because the magnetisation is allowed to fully recover (i.e., fully recharge) between each excitation of the magnetisation. While this seems good, what do the images produced with these TRs look like? What TR choice produces better WM/GM contrast? How much time does such a scan take?

*Images produced with very long TRs and short TEs do not contain any T1 (or T2) weighting. If the magnetisation in all tissue types is allowed to fully recover, then the ability to differentiate tissues based on the rate of recovery is lost. Therefore, these images show only small differences in signal strength across different tissue types (only reflecting proton density), and are more or less uniformly bright looking. These types of images only tell you how much water there is at each spatial location.*

*Better T1-weighting can be produced by reducing the TR so that full  $M_z$  recovery of*

*all tissue types is not achieved. In a sense, you want to measure your signal before the slow recovering (long T1) tissues have time to “catch up”. Reasonable TR values range from 20 – 1000 ms (ignoring SNR and scan duration considerations), with scan times ranging from 8 mins to ~7 hours! Good CNR has to be balanced with appropriate scan times.*

**Question 1.2\***

Given a TR of 100ms, a square field-of-view (192x192) and assuming a 2DFT (single line) readout with 128 slices, estimate the total imaging time required.

*ScanTime =  $N_{PE} \cdot TR \cdot N_{SLICES}$  giving 2458 seconds, or about 41 mins.*

Does this time seem reasonable? If not, what might you change?

*You could reduce the TR, or the field of view (which would lead to larger voxels).*

As you learned in last week’s tutorial, the steady-state signal (i.e., the signal available after long series of excitations) is often maximised at a flip angle  $< 90^\circ$  when the TR is not long relative to T1. Recall that the Ernst angle  $\theta_E$  is defined:

$$\cos(\theta_E) = e^{-\frac{TR}{T1}}$$

In addition to T1-contrast levels and scan time, signal-to-noise ratio (SNR) and contrast-to-noise ratio (CNR) must also be considered when selecting imaging parameters.

*Fix TR = 20 ms, and demonstrate the effect of changing flip angle on the image SNR and CNR. Discuss the results in the context of changing steady-state signal levels that arise from tissues with different T1 values.*

*Demonstrate the use of `simSignalvFlip` and `simContrastvFlip` to explore these effects.*

```
>> T1 = 600;  
>> simSignalvFlip;  
>> T1 = 900;  
>> simSignalvFlip;  
>> T1pair = [600, 900];  
>> simContrastvFlip;
```

**Question 1.3**

What flip angle maximises the contrast between tissues with T1 values of 600 ms and 900 ms when a TR of 20 ms is being used?

*This value can be read off the plot generated by `simContrastvFlip` (students may have to zoom in). A flip angle of  $\sim 23^\circ$  will maximize WM (600 ms T1)/GM (900 ms T1) signal contrast. Note that this contrast maximum is distinct from the relative signal maxima associated with WM alone ( $\sim 15^\circ$ ) and GM alone ( $\sim 12^\circ$ ).*

**Question 1.4\*\***

How are noise levels affected when the acquisition flip angle is modified?

*They aren't! However, other imaging parameters can affect noise, for example increasing bandwidth would increase noise (but reduce the acquisition time, as larger bandwidth means applying stronger gradients).*

---

## Part 2 – T2 Contrast [20–25 minutes]

*You may want to instruct students to close all their open figure windows to reduce clutter*

```
>> close all
```

T2-weighting is conceptually identical to T1-weighting, except now the objective is to maximise signal differences based on differences in tissue T2 values (sometimes also seen as  $T_2$ ), instead of T1 values. Recall that T2 is distinct from T1 in that it governs the decay or loss of transverse magnetisation (as opposed to the recovery of longitudinal magnetisation). The T2 value indicates the speed of this decay, where tissues with short T2 decay quickly, and tissues with long T2 stay around longer and take longer to decay away.

*This can be easily demonstrated by opening another simContrast window and exploring the lower lefthand subplot.*

#### Extra Information

For the mathematically inclined, relaxation of the transverse magnetisation is governed by the following equation:  $M_{xy}(t) = M_z \sin \theta e^{-t/T_2}$ , where  $M_z$  is the longitudinal magnetisation at  $t = 0$  and  $\theta$  is the flip angle.

An important point when trying to optimise image contrast is that T1 and T2 weighting are *always* present to varying degrees in acquired images. In the previous section, TEs were set to very low values (or zero) to minimize T2-weighting effects. Similarly, when trying to maximise T2 weighting it is beneficial to minimise the effects of T1 differences in tissue. This can be done in 2 ways:

- Increasing TR
- Decreasing flip angle

*Explore the combined effects of increasing TR and decreasing the flip angle on T1 contrast. Find a combination that leaves the image with little T1-weighting, and high SNR. Something like TR = 200, flipAngle = 12 works well. Keep TE = 0 for this. (Note that in practice, flip angles of around 90° are more commonly used, to allow for the acquisition of multiple lines in k-space following each excitation.)*

In the same way that changing TR changed T1-contrast in the signals, changing TE modifies the amount of T2-contrast in the resulting image.

*Demonstrate the change in contrast with TE.*

```
>> TE = 0;
>> simContrast
>> TE = 100;
>> simContrast
```

*Note that although the SNR drops, the CNR increases at the larger echo time.*

*Use simContrastvTE to demonstrate the variability in image contrast with different choices of TE.*

```
>> T2pair = [80, 100];
>> simContrastvTE
```

#### Question 2.1

At approximately what TE is the contrast between WM and GM at 1.5 T maximised? What implications does this have for minimum scan times for a simple 2DFT (single line) sequence?

*Signal contrast is maximised when TE ~ 90 ms. More generally, given two tissue types with different T2 values, contrast will be maximised when TE is somewhere in between the two T2s. Because TE < TR (in most cases), imaging at this TE using standard methods would lead to extremely long scan times.*

#### Question 2.2\*\*

In general, given two tissue types with different T2 values, derive an expression for the TE that gives maximum contrast, assuming TR >> T1.

*Let  $S_A = S_{0,A} e^{-TE/T_{2A}}$  and  $S_B = S_{0,B} e^{-TE/T_{2B}}$ , then by solving  $\frac{d(S_A - S_B)}{dTE} = 0$ , you find that*

$$TE_{OPT} = \frac{\log\left(\frac{T2_B \cdot S_{0,A}}{T2_A \cdot S_{0,B}}\right)}{\left(\frac{1}{T2_A} - \frac{1}{T2_B}\right)}$$

### Discussion 2.3\*\*

Often, SNR is said to be proportional to the square root of the total acquisition time ( $T_{acq}$ ). Consider two pulse sequences that are identical in all respects, except one has double the TR of the second. If we were to apply the acquisition time “rule”, the SNR of the second sequence should be  $\sqrt{2}$  higher than the first. Assuming a homogenous sample (single  $T1$ , single  $T2$ ), is this relationship valid? What happens to the SNR when the TE is varied to be as short as possible within a given TR, or as long as possible?

- *The relationship is not valid as SNR is actually proportional to  $\sqrt{\text{readout time}}$ , not total acquisition time. Be careful here, as changes to TR will affect the Ernst angle and steady-state signal levels.*
- *Similarly here, moving TE doesn't change readout duration, but it does affect the signal level depending on T2.*

### Extra Information

Remember from the lecture that  $T2$  relaxation is an intrinsic process that occurs in any sample, regardless of how perfect (or imperfect) your instrumentation or experimental setup is. In practice, static magnetic field inhomogeneity across a voxel can cause an apparent enhancement of  $T2$  decay or transverse signal loss. We call  $T2^*$  (pronounced “tee two star”) the parameter that takes into account both intrinsic  $T2$  and field inhomogeneity effects, where  $T2^* < T2$  (i.e., inhomogeneity always causes faster decay).

## Part 3 – Inversion Recovery [20 minutes]

Inversion pulses are just like excitation pulses, except their purpose is to flip the longitudinal ( $M_z$ ) magnetisation upside down along the z-axis. Once this occurs, the magnetisation immediately begins to grow back to its positive equilibrium value, governed by its  $T1$  parameter. *Preparing* the magnetisation this way allows more flexibility in dealing with tissues with different  $T1$ s, including the ability to selectively “null” tissues. The sequence parameter TI controls the timing between the inversion pulse and the sequence excitation pulse.

Approximate Relaxation Times at 3.0T		
	T1 (ms)	T2 (ms)
White Matter (WM)	830	80
Grey Matter (GM)	1330	110
Cerebrospinal Fluid (CSF)	4000	2000

Table 2

*Notice the differences in T1 and T2 times between 1.5 and 3T.*

Open a new `simContrast` figure with the following parameters:

```
>> fieldStrength=3;
>> TR = 2000;
```

```
>> TI = 860;
>> TE = 5;
>> flipAngle = 90;
>> useInversion = 1;
>> simContrast
```

*Take some care in going over the upper left subplot, particularly the relationship between the excitation and inversion. Also note that the above parameters were chosen to null CSF.*

### Question 3.1

Assuming  $TR = 1000$  ms,  $TE = 5$  ms, and flip angle  $= 90^\circ$ , find the approximate  $TI$  that corresponds to nulling of white matter signals.


*A  $TI \sim 312$  ms corresponds to a white matter null.*

### Question 3.2

Assuming  $TR = 1000$  ms,  $TE = 5$  ms, and flip angle  $= 90^\circ$ , find the  $TI$ 's (there are two) that correspond to WM and GM having identical signal intensities.

*$TIs$  of 192 and 335 ms lead to identical WM/GM signal magnitudes. In the former case, this is where the  $M_z$  recovery curves for both WM and GM intersect, and in the latter case, it is where the WM and GM curves have equal magnitude but opposite polarity (WM positive, GM negative).*

### Question 3.3\*\*



Find an expression for null inversion time of a tissue with an arbitrary  $T_1$ , using a general steady state sequence with arbitrary  $TR$  and flip angle.

Hint: the longitudinal magnetization  $M_\theta$  after a flip angle  $\theta$  is given by

$$M_\theta = M_0 + (M_z \cos \theta - M_0)e^{-t/T_1}$$

where  $M_z$  is the longitudinal magnetization immediately before applying the flip angle and  $M_0$  is the equilibrium magnetization.

First find the null inversion time for an arbitrary tissue with an arbitrary starting longitudinal magnetisation  $M_z$ .

Using the above equation set  $M_\theta=0$ ,  $\theta = 180$  degrees and solve for  $TI$ :

$$M_0 - (M_0 + M_z)e^{-TI/T_1} = 0$$

$$TI = T_1 \cdot \log \left( 1 + \frac{M_z}{M_0} \right)$$

Secondly, for a general steady state sequence with an inversion preparation (as in the figure above) compute the steady state longitudinal magnetisation  $M_z$  by setting  $M^n = M^{n+1}$  (as in tutorial 2, question 2.4):

$$\frac{M_z}{M_0} = \frac{(1 - e^{-(TR-TI)/T_1} + \cos(\theta) e^{-(TR-TI)/T_1} - \cos(\theta) e^{-TR/T_1})}{1 + \cos(\theta) e^{-TR/T_1}}$$

Since the magnetisation for the tissue we are considering is defined to be 0 after time  $TI$ , the expression for  $\frac{M_z}{M_0}$  becomes a simple  $T_1$  recovery:

$$\frac{M_z}{M_0} = 1 - e^{-(TR-TI)/T_1}$$

The final expression relating  $TI$ ,  $TR$ , and  $T_1$  is

$$TI = T_1 \cdot \log \left( \frac{2}{1 + e^{-TR/T_1}} \right)$$